



PC TABLETS FOR MICROSCOPY



The new **MICROPAD-7** and **MICROPAD-10** are compact and comfortable PC Tablets with 7" and 10" displays and a built-in 5 MP camera that can be directly connected to the microscope using the step C adapter to the third output of the microscope. Multiple features make them excellent for scientific teaching.

MICROPAD-7




MICROPAD-10





Excellent imaging solution provider



TECHNICAL SPECIFICATIONS

Cod.	MICROPAD-7	MICROPAD-10
Display	LCD 7" TOUCH SCREEN backlit	LED 10" TOUCH SCREEN Multitouch
Display resolution	1024x600	1024x768
Operating system	ANDROID 2.3	WINDOWS 7
CPU	A8 800Mhz	Intel ATOM N570 (1.66GHz)
Memory	RAM 512MB DDR2-400	2GB 1066MHz DDR3 SDRAM
Graphic Card	UMA Integrated 200MHz/DirectX9	UMA Integrated 200MHz/DirectX9
Bluetooth	Supported	Supported
Outputs	Mini USB 2.0 Micro HDMI Micro SD SIM card 3G External	USB 2.0 (X2) HDMI Micro SD SIM card 3G External External interface LAN & VGA
Wi-Fi	Wireless data network, WLAN 802.11b/g	Wireless data network, WLAN 802.11b/g
Ethernet Connection	/	Onboard 10/100 Mbps
Monitor Inclination	Up to 90°	Up to 90°
Working temperature	-10 ° to 35 ° C	-10 ° to 35 ° C
Power Supply	DC5V 250mA	DC100-240V AC 50/60Hz 19V/40W P.A.
Battery	3300mAh Rechargeable - Duration 4 HOURS	3300mAh Rechargeable - Duration 4 HOURS
Built-in Camera for Microscope	5 Megapixel	5 Megapixel
Sensor	CMOS 1/2.5"	CMOS 1/2.5"
Image Previewing	30 FPS full resolution	30 FPS full resolution
Video Format	480P/15FPS	480P/15FPS
Selectable resolution	2592x1944/1280x960/1024x768/640x640	2592x1944/1280x960/1024x768/640x640
Lens mount	C-mount	C-mount
Software	LIVE	LIVE - MEASUREMENTS - REPORT



Pulse Healthcare Technology
House: 11/1, Shahid Minar Road
Kallyanpur, Dhaka-1216, Bangladesh.
Mobile: +8801708008061
e-mail: info@pht.com.bd
web: www.pht.com.bd



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

